




## PleuraFlow® Active Clearance System® In Use



**The PleuraFlow Active Clearance System® is not MRI compatible.  
Do not take patient to MRI with PleuraFlow System® in place.**




For Clinical Support call 1-844-257-3569 ext2

DC101-C



## PleuraFlow® Active Clearance System® In Use



**The PleuraFlow Active Clearance System® is not MRI compatible.  
Do not take patient to MRI with PleuraFlow System® in place.**



For Clinical Support call 1-844-257-3569 ext2

DC101-C